

(19) World Intellectual Property
Organization
International Bureau



(43) International Publication Date
17 June 2004 (17.06.2004)


PCT

(10) International Publication Number
WO 2004/049841 A3

- (51) International Patent Classification⁷: **A61F 5/03** (74) Agent: **BARKER, Brettell**; 12-10 Priest's Bridge, London SW15 5JE (GB).
- (21) International Application Number: PCT/IB2003/006395 (81) Designated States (*national*): AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BW, BY, BZ, CA, CH, CN, CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, EG, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NI, NO, NZ, OM, PG, PH, PL, PT, RO, RU, SC, SD, SE, SG, SK, SL, SY, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, YU, ZA, ZM, ZW.
- (22) International Filing Date: 27 November 2003 (27.11.2003) (84) Designated States (*regional*): ARIPO patent (BW, GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZM, ZW), Eurasian patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European patent (AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HU, IE, IT, LU, MC, NL, PT, RO, SE, SI, SK, TR), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG).
- (25) Filing Language: English
- (26) Publication Language: English
- (30) Priority Data: 0227838.0 29 November 2002 (29.11.2002) GB
- (71) Applicant (*for all designated States except US*): **QUALITEAM S.A.S.** [IT/IT]; Casale Nasso Sopra 15, I-10010 Chiaverano (IT).
- (72) Inventor; and
- (75) Inventor/Applicant (*for US only*): **HJORTH HANSEN, Doris** [DK/IT]; Casale Nasso Sopra 15, I-10010 Chiaverano (IT).
- Published:
— with international search report

[Continued on next page]

(54) Title: POST-OPERATIVE VEST



(57) **Abstract:** The invention relates to a variable tension chest support for use by patients having undergone thoracic surgery. The chest support may be worn directly on the skin of a patient and comprises a band of stretchable material (1) and one or more grips (6) which are wholly or partially concealed within pockets (4) with which the patient may alter the tension of the chest encircling band. The device allows a patient to modify the level of tension applied to the chest, thus managing their recovery process, while providing a basic level of support at all times the support is in use.